NO Bruksanvisning	3 - 6
EN User manual	7 - 10
DE Bedienungsanleitung	11-14
FR Manuel d'utilisation	15-18
NL Gebruiksaanwijzing	19-22

A 814728

M✓ S✓ X

B 814059

M✓ S✓ X

C 814025

M✓ S✓ X✓

D 814765

M✓ S✓ X✓

E 814622

M✓ S✓ X

F 814046

M✓ S✓ X

G 814009

M✓ S✓ X

H 814043

M✓ S✓ X✓

I 814024

M✓ S✓ X✓

J 814023

M✓ S✓ X✓

K 814141

M✓ S✓ X✓

L 814014

M✓ S✓ X✓

M 814032

M✓ S✓ X✓

N 814729

M✓ S✓

814730


X✓


TOPRO TROJA M
TOPRO TROJA S
TOPRO TROJA X

814749
814742
814741


1


2


3


4


5


6

7

8

9

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Welcome as a user of a TOPRO PRODUCT

Every TOPRO product is developed, designed and manufactured in Norway. Due to advanced manufacturing by using the best equipment we are able to offer great accuracy and maximum quality in all levels of our production. This product conforms to the Medical Device Directive 93/42/EEC. It is tested and approved according to EN-ISO 11199-3:2005.

On receipt of the goods, please check your product, if you have any queries please contact your mobility shop or dealer immediately. The same applies if you have general questions about the product. If you have problems reading this user manual an electronic version is available on www.topro.co.uk

Guarantee // Expected life span of the product

TOPRO TROJA WALKER is guaranteed free of faults and defects for a period of 7 years.

An exception is if the damage is caused by incorrect use or parts that are exposed to natural wear and tear (such as e.g. brake block, brake cable, wheel, seat, basket and handgrip). For repairs during the guarantee period, please contact your local mobility shop or dealer. The guarantee is void if unauthorised spare parts and accessories have been or are being used on the product.

The expected life span of the product is estimated to be 10 years, provided that it is used in a regular way according to this user manual and its safety and maintenance instructions.

The red numbered squares below refer to the numbered illustrations on the inside of the cover.

Get to know your TOPRO PRODUCT **1**

- | | | |
|-------------------------------|------------------------------------|-----------------------------------|
| A Forearm support | F Seat bars | K Basket |
| B Brake handle | G Seat | L Side frame |
| C Brake lever | H Lock clamp | M Brake block / Edge guard |
| D Forearm support tube | I Strap for folding | N Tilt function |
| E Brake wire | J Lock wheel for the handle | O Supporting cord |

P The label on the rear side of the tube under the seat shows the product name/model and

- Maximum weight of user
- Maximum width of walker
- Maximum height of walker

A separate label on the tube indicates the serial number (SN) and production year and month.

Measurements

Model/size	TOPRO TROJA M	TOPRO TROJA S	TOPRO TROJA X
Colour	Dark grey	Silver grey	Silver grey
Recommended user height	150–185 cm	130–170 cm	130–160 cm
Maximum weight of user	135 kg	125 kg	80 kg
Maximum load basket	5 kg	5 kg	5 kg
Height of support points	980–1170 mm	840–1000 mm	840–960 mm
Maximum length	650 mm	650 mm	650 mm
Maximum width	630 mm	630 mm	530 mm
Between support points	450 mm	450 mm	350 mm
Seat height	620 mm	540 mm	500 mm
Turning diameter	840 mm	840 mm	790 mm
Diameter wheels	200 mm	200 mm	200 mm
Width wheels	36 mm	36 mm	36 mm
Diameter hand grip	35 mm	35 mm	35 mm
Grip distance driving brake	65 mm	65 mm	65 mm
Height when folded	1020 mm	880 mm	880 mm
Width when folded	330 mm	330 mm	330 mm
Length when folded	650 mm	650 mm	650 mm
Weight walker	10,1 kg	9,8 kg	9,3 kg

Intention

The TOPRO TROJA WALKER shall support users with poor balance and/or poor mobility. The integrated seat provides an opportunity for the user to rest if needed, and the fabric basket is capable of carrying items up to a total of 5 kg. The product is not suitable for persons with low strength in arms, with very poor balance, or with considerable cognitive disabilities. The walker is designed and approved to be used indoors on plain foundation. It is developed to be pushed, not pulled.

Important to know about safety before use

- Take care when folding and unfolding, do not trap your fingers.
- Check the brakes every time before use and that the walker locks into position when opened.
- Not to be used by persons over maximum user weight (see table).
- Do not use to transport heavy loads or persons.
- The walker is not a toy.
- Make sure that the lock wheel for height adjustment for the handles is locked tight.
- Be aware! Some parts of the walker might feel hot or cold when exposed to extreme temperatures.
- Do not alter the product as this might put your safety at risk and the guarantee will not be valid.

Important to know about safety when in use **2**

- When stationary or parked, always use the parking brakes on both sides **A**.
- Before sitting on the seat, make sure the parking brakes on both sides are in the locked position **A**.
- Do not lean or bend sideways or over the back of the walker whilst in the seated position. Do not tilt the walker, to avoid risk of falling. Due to the same reason, bags etc. must never be hung on the handlebars.
- The driving brakes can be held in a little **B** to reduce speed, but must not be in locked position when the walker is in motion.
- To give the right support and a safety, make sure the walker is kept as close to your body as possible **C**.
- Walk straight, look forward and hold your backside in **D**. That is the easiest way to keep your balance.
- Do not push the walker too far from your body as this can cause unnecessary strain and you could fall.
- Do not use the walker on stairs. Be careful when descending thresholds etc.
- Be careful of flared skirts and trousers as they could get caught between the brake block and wheel, which could lead to falling. We therefore recommend not wearing flared/wide clothing.
- For safety reason it is prohibited to sit on the rollator during public transportation. Place your walker securely and use the parking brake.

Assembling / adjusting

TOPRO TROJA is delivered assembled. Before use, you must unfold the walker. The height of the forearm support tubes will then need adjusting to suit the user's needs.

Unfolding / folding **3**

Unfolding: Loosen the strap between the forearm supports **A**. Loosen the lock clamp that is over the seat bars **B**. Press the seat bars down and out **C**. The side frames will slide apart. Press down the seat bars. You will hear a clicking sound. The walker is now locked in its open position. Fasten the strap **A**.

Folding: Pull the strap on the seat **D** and press the side frames together. Lock the clamp **E** over the seat bars.

Adjusting the handles and height adjustment of the forearm supports **4**

Adjusting the handles: Unscrew the lock wheel **A** until the handle can be moved, adjust to the wanted position and tighten the lock wheel again. It is done separately on each side. It is not possible to pull the handles all the way out; there are physical stops both ways. Height adjustment: Unscrew the lock wheel **B** (approx. 4-5 turns). Adjust the forearm support tubes to a suitable height, then pull them slightly up and down **C** until you hear a clicking sound. Tighten the lock **B** again. Adjust both handlebars to the same height.

Recommended height: Stand with relaxed shoulders with your elbows in a 90 degree angle. **2**

Height adjustment using the memory function **5**

You have the opportunity to set the height of the forearm supports at a standard setting fitting to your height so that you easily can find the same height again. This can be necessary when the walker is being transported or stored as compact as possible. Decide the height of the forearm support tube. Count visible holes at the back of the tubes. Unscrew the lock wheel **A**. Pull the tube **B** out of the frame. The attached part **C** will now follow the tube. Move the attached part **C** on the tube **B** so there are equal holes as counted earlier. Fit the memory chip **D** on the tube in the first available hole under the attached part **C**. Push the tube **B** down and into the side frame **E** until the attached part is locked. Screw the lock wheel **A** half in. Pull the tube **B** up until it stops, and then tighten the lock wheel **A** properly. From now on, when the forearm support tubes are in the lowest position, they will stop in the pre-set position automatically when pulled up.

Adjusting / checking and maintaining the brakes **6**

Adjusting / checking: Loosen the nut **A** beside the brake. Tighten the brake wire **B** with the adjustment screw **C**. Adjust both brake blocks so the distance to the wheel surface is approximately 1 mm **D**. Tighten the nut **A**. Make sure the brakes are not adjusted too tight and pressing into the rubber when the parking brake is on.
Maintaining: Clean the brake block **E** and the brake spring **F** frequently. The wheels can be changed if they are worn out or damaged.

Using the brakes **2**

Make sure the brakes work before each trip. Each brake handle **A** works on each of the rear wheels. The walker wheels are made of a special soft rubber, which gives an optimal braking capacity on slippery floor. Wrong use can cause unsafe situations and damage to the walker: **The walker must not be driven or pushed with force when the parking brakes are on or when the wheels are partly locked using the brakes. This indicates incorrect use and can cause wear and tear which makes the tires and brake quality weaker.**
Parking brakes: Push both brake levers **A** forwards until they lock. Both rear wheels are now locked. Pull the brake levers back to release the parking brakes. *Driving brakes:* Squeeze both brake handles **B** upwards to reduce speed. **Do not use continuously while the walker is in motion.**

Seat and basket **7**

The parking brakes must be on before using the seat **A**. Push both brake levers **B** forwards until they lock. Both rear wheels are now locked. Sit with your back towards the basket **C**. The basket **D** can carry goods up to 5 kg. It can be removed and washed separately and gently (hand wash 40°C).

Edge guard and tilt function **8**

Edge guard: The brake blocks **A**, which are in front of the rear wheels, prevent the rear wheels from catching on corners, door frames etc. **Tilt function:** To use the tilt function use the front of the foot and press down on the tilt pedal **B** which is beside the rear wheel. This helps to tilt the walker over thresholds, edges etc.

Carrying the walker **9**

The walker must not be carried by the seat strap or the brake wires. Fold the walker together as shown on the illustration. Carry the walker by the folded seat bars **A**.

Transport

The walker should be transported in the folded and locked position (see paragraph *Unfolding / Folding*). Ensure that the walker cannot unfold unintentionally. Be careful when loading the walker into a vehicle, giving attention to securing of the item.

Storage

The walker should be stored in an upright position. It can be stored in and outside, if stored outside then under cover.

Material // Recycling

The walker is made of plastic coated aluminium pipes, plastic couplings and material of plastic and polyester. Most parts can be recycled. Dispose the walker according to the applicable regulations in your country. For information, please contact your public authorities.

Cleaning

The walker can be steamed, hosed and washed with normal household detergents. Do not use scrub detergents, steam- or high-pressure cleaners.

Disinfection

Disinfection must only be carried out by authorized personnel and by wearing sufficient protective equipment. The surface must be cleaned with a disinfection that contains 70-80% ethanol. We advise against disinfection that contains chlorine and phenol. *The manufacturer cannot be held responsible for any damage or wounds that might occur by the use of harmful disinfection or disinfection carried out by unauthorized personnel.*

Inspection // Maintenance // Re-use

It is recommended that general inspection and maintenance is carried out and that the items below are checked on a regular basis (how frequent is determined by how and how often the walker is used): Tubes, frame and screws, forearm supports and handles, brakes, wheels, seat and accessories. This also applies when the walker is to be prepared for re-use. Please see detailed instructions in this user manual about maintenance of the brakes and cleaning and disinfection of the walker.

If the walker does not work properly

Make sure the walker is unfolded correctly and that both brake handles are in expected position.

Accessories

For pictures and TOPRO article numbers see the inside of the cover. Accessories mounted on the walker may influence stability, we therefore advise you to use them carefully.

- | | |
|---|--|
| A Tray | H Bottle holder |
| B Anti slip mat for tray | I Nameplate |
| C Crutch holder | J Lock |
| D Back support | K Flower stickers, set with 4 stickers |
| E Shopping bag | L Guiding handle for e.g. visually impaired |
| F Rear bag with zipper | M Drag brake (pair) |
| G Basket for oxygen bottle (only compressed, not liquid) | N Stabilizer bar |

The accessories are to be ordered separately to give the possibility to compose a TOPRO TROJA WALKER that suits the individual users need. Contact your mobility shop, dealer or TOPRO for an updated overview of the accessories, or visit our homepage www.topro.co.uk

In case of faults or questions please contact your mobility shop or dealer.

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