

SEAT SYSTEMS For REAL chairs







ERGOMEDIC – ADVANCED ERGONOMICS

With various backs, seats, castors, lifting columns and accessories, it is easy to put a functional, comfortable chair together. To further increase sitting comfort, Mercado Medic has collaborated with registered physiotherapist, Bengt Engström, to develop the ErgoMedic seat system for the REAL 9000. ErgoMedic consists of three interacting parts: Seat, Back and the Armrest.

Backrest

The ERGOMEDIC back rest and seat work together to provide the most comfortable, and effective sitting position. The shape and angle of the high back has been improved to create more space for the pelvis and its softer parts. The lower edge of the back panel now has a gentle 'U' shape and is softer. This, combined with the shape of the seat, provides better sitting depth and a more comfortable position for the lower back and pelvis.

Seat

For optimal sitting comfort, the ERGOMEDIC® seats have the right form and softness to ensure optimal stability and comfort, higher at the thigh than the buttocks. Its form distributes pressure evenly and effectively, and prevents the user from sliding forward when the seat is angled. The back edge of the seat is set low to increase the sitting depth.

This means that the seat bone nodules can move backwards when the users stretch out forwards. When the pelvis is positioned at the very back of the seat, and well supported by the backrest, the users legs are also stabilised.

Coxit seat

Mercado's Coxit seat has the same high ergonomic quality as the ERGOMEDIC seat. The Coxit seat is divided with individually adjustable flaps. It is used primarily for the alleviation and rehabilitation of orthopaedic injuries to hips or thighs.

Armrest

The armrest is an extremely important part of the ERGO-MEDIC system, so we have also developed these to be more ergonomic and comfortable. The edges have been bevelled, both on the inside and the outside to minimise the pressure against the underarms during prolonged contact.

The armrest panels are hollowed out at the back in order to stabilise the elbows. At the front they have been designed to provide a good grip when standing up and transferring pressure.

The armrests are covered with soft imitation leather that is more pleasant to the touch and is less slippery if the armrest becomes damp.



Bac	krest ErgoMedic, width x depth	Article No.
1	Backrest low ErgoMedic 27x22 cm	210466
2	Backrest low ErgoMedic 30x26 cm	210464
3	Backrest low ErgoMedic 33x26 cm	210462
4	Backrest low ErgoMedic 38x26 cm	210460
5	Backrest high ErgoMedic 27x32 cm	210467
6	Backrest high ErgoMedic 30x34 cm	210465
7	Backrest high ErgoMedic 33x36 cm	210463
8	Backrest high ErgoMedic 40x34 cm	210461
9	Backrest high ErgoMedic 42x43 cm	210468
Sea	t ErgoMedic, width x depth	Article No.
10	Seat ErgoMedic 29x32 cm	150468
10	Seat ErgoMedic 33x35 cm	150466
10	Seat ErgoMedic 36x39 cm	150465
10	Seat ErgoMedic 36x43 cm	150467
10	Seat ErgoMedic 40x39 cm	150464
10	Seat ErgoMedic 40x43 cm	150463
10	Seat ErgoMedic 40x48 cm	150469
10	Seat ErgoMedic 44x43 cm	150462
10	Seat ErgoMedic 44x48 cm	150460
10	Seat ErgoMedic 48x48 cm	150461
Sea	t Coxit ErgoMedic, width x depth	Article No.
11	Seat Coxit Int ErgoMedic 40x43 cm	160463
11	Seat Coxit Int ErgoMedic 40x48 cm	160469
11	Seat Coxit Int ErgoMedic 44x43 cm	160462
11	Seat Coxit Int ErgoMedic 44x48 cm	160460
11	Seat Coxit Int ErgoMedic 48x48 cm	160461

EXPLANATION OF THE STRUCTURE OF THE ARTICLE NUMBER

Article No: 150460. "ErgoMedic seat, plush grey 44x48 cm"

Digit 1 = Type of upholstered article Digits 2-3 = Seat system Digit 4 = Material Digit 5 = Colour Digit 6 = Size



ERGOMEDIC[®] PLUS – ERGONOMIC DESIGN THAT IS OUT OF THE ORDINARY

Backrest

The ERGOMEDIC® PLUS back is available in several variants. The higher backrest provides more shape around the pelvis and thoracic spine, whereas the lower backrest has the same shape, but supports the pelvis and the lumbar spine.

The lower part of the high backrest is shaped so that you can sit with pelvic support even when tilted backwards. At the top in the centre, between the shoulder blades, a soft extra support has been integrated. This to enable the user to sit upright and stable against the backrest, whilst at the same time leaving enough room for free movement of the shoulders.

The high backrest is also available with extra side wedges. The low backrest has a similar shape to the high one, but is designed for supporting the pelvis and lumbar spine.

Seat

The biggest requirement of the seat is to be stable and comfortable. ERGOMEDIC® PLUS seats have the right form and softness to ensure optimal stability and comfort, higher under the thigh than the buttocks. Its form distributes pressure evenly and effectively, and prevents the user from sliding forward when the seat is angled. When the user sits at the back of the seat, so that the pelvis is firmly placed against the lower parts of the back support, their legs are also stabilized. The higher rounded shape on either side of the thighs means that both the pelvis and the legs are stabilized.

Armrest

The user should be able to switch between active, and relaxed positions. The ERGOMEDIC® PLUS armrests are shaped so the arms are supported from the sides, which means that the upper body is more relaxed. They are also raised beneath the wrists to provide extra relief.

At the front, the armrest has a rounded convex shape, which makes it comfortable to rest the hand on , while providing a good grip when you want to get in or out of the chair

Neckrest

The neckrest provides support primarily against the rear, lower part of the head. The lower part of the neck rest is built up with soft padding on

each side, providing support for the lower skull and neck. The upper part is flatter, which leaves room for the user to turn their head. The shape and properties of the neck rest means that the head is well stabilised, even when the user is sitting upright.



Article No. 980825 251460

251465

251463

251468

251469

252465

252463

252468

252469

254463

254468

Dackrest Ergoweald PL05, width x depth		
1	Backrest ErgoMedic Plus	35x20 cm
1	Backrest ErgoMedic Plus	40x20 cm

Poolstost ErgoModio DLUS width x donth

- Backrest ErgoMedic Plus 30x34 cm 2
- 3
- Backrest ErgoMedic Plus 35x43 cm
- 4 Backrest ErgoMedic Plus 42x43 cm
- 4 Backrest ErgoMedic Plus 46x43 cm
- 5 Backrest EM Plus with wedge 30x34 cm 6 Backrest EM Plus with wedge 35x43 cm
- 7 Backrest EM Plus with wedge 42x43 cm
- 7 Backrest EM Plus with wedge 46x43 cm
- Backrest EM Plus Comfort 35x43 cm 8
- 8 Backrest EM Plus Comfort 42x43 cm
- 8 Backrest EM Plus Comfort 46x43 cm
- 254469 8 Backrest EM Plus Comfort w wedge 35x43 255463
- 8 Backrest EM Plus Comfort w wedge 42x43 255468
- 8 Backrest EM Plus Comfort w wedge 46x43 255469

Seat	t ErgoMedic PLUS, wi	dth x depth	Article No.
9	Seat ErgoMedic Plus	29x32 cm	151468
9	Seat ErgoMedic Plus	33x35 cm	151466
9	Seat ErgoMedic Plus	36x39 cm	151465
9	Seat ErgoMedic Plus	36x43 cm	151467
9	Seat ErgoMedic Plus	40x39 cm	151464
9	Seat ErgoMedic Plus	40x43 cm	151463
9	Seat ErgoMedic Plus	40x48 cm	151469
9	Seat ErgoMedic Plus	44x43 cm	151462
9	Seat ErgoMedic Plus	44x48 cm	151460
9	Seat ErgoMedic Plus	48x48 cm	151461

Sea	t Coxit ErgoMedic PLUS, width x depth	Article No.
10	Seat Coxit Int EM Plus 40x39 cm	161464
10	Seat Coxit Int EM Plus 40x43 cm	161463
10	Seat Coxit Int EM Plus 40x48 cm	161469
10	Seat Coxit Int EM Plus 44x43 cm	161462
10	Seat Coxit Int EM Plus 44x48 cm	161460
10	Seat Coxit Int EM Plus 48x48 cm	161461
Nec	krest, width x height	Article No.

11 Nekrest cushion EM Plus 30x23 cm

Armrest

- 12 Armrest ErgoMedic Plus Right 42 cm
- 12 Armrest ErgoMedic Plus Left 42 cm
- Armrest ErgoMedic Plus Right 30 cm 12
- Armrest ErgoMedic Plus Left 30 cm 12

★ The ErgoMedic Plus Comfort back is designed for maximal comfort. It has an adjustable (pump) support for the inward curve of the lower part of the back.

270465

- Article No. 243380 243385
- 243381 243386



OTHER UPHOLSTERED PARTS

Trunk support

The trunk support is used to produce more upright, active sitting. Mercado has three types: Curved and adjustable, straight, and curved trunk support.

Neck rest

The neckrest provides support for the head in order to achieve the best sitting position and to provide relief for the neck and spine even for users with poor body control. It can easily be adjusted to achieve the best setting for the individual user. Available in flat, curved or ERGOMEDIC® PLUS versions.

Cover

The covers protect the upholstery from dirt and dust, as well as extending the upholstery's lifetime. Covers are available in different colours for a variety of combinations.

Incontinency covers are a quick and easy solution to protect seats against moisture.

Armrest

The Hemiplegi armrest is suitable for partially paralysed, or users with restricted control. The armrest panel has been hollowed out, which makes it easier to keep the paralysed arm on the rest. The armrest panel also rotates, which allows for the best possible position. The mobile armrest panel is upholstered with extra foam for optimised comfort.

Side support

The side supports help the user maintain stability and position. These can also be used as thigh supports, so that the users feet are better supported on the footrests or the floor, and prevents their legs from sliding apart. It can be adjusted for height and width, or turned backwards and forwards. Available in three sizes: S, M and L.

Slow Recovery seat and others

Slow Recovery seats are manufactured with two materials in the cushioning, which provides optimum support and comfort to the user. The upper layer comprises of a memory foam that moulds itself, and provides good pressure distribution and increased sitting comfort, while the lower layer raises the seat bone nodules from the wood plate. Even if the user does not sit in the chair for a sufficiently long period of time, the seat cushion adjusts itself so that a cavity is formed where the person sits. This material is available for all upholstered parts.

The wooden seat disc should be used when casting or retrofitting of its shell. The seat with anti-slip is recommended when using individually adjusted seats. Seat with hygiene fabric should be used by users with incontinence problems. To be completed with a cover in plush.



Tru	nk support, width x heig		Article No.	Neckrest, width x	-	Article No.
1	Trunk support cushion		280460		ion Straight 29x16 cm	
2	Trunk support curved a	nd adjustable	801775	4 Neckrest cush	ion Curved 31x18 cm	271460
Arm	nrest		Article No.	Article No. for sea	t angle adjustment	
6	Armrest ErgoMedic 25	cm	240382	240387		
6	Armrest ErgoMedic 30	cm	240381	240386		
6	Armrest ErgoMedic 35	cm	240380	240385		
6	Armrest ErgoMedic 40	cm	240383	240388		
7	Armrest Mobil 25 cm		245382	245387		
7	Armrest Mobil 30 cm		245381	245386		
7	Armrest Mobil 35 cm		245380	245385		
7	Armrest Mobil 40 cm		245383	245388		
8	Armrest Hemiplegi 45x	12 cm	244382	244387		
9	Wooden seat disc		15090x	Covers		
10	Seat with anti-slip		15061x		patterned fabric and in n to www.mercado-me	
Side	e support, width x heigl	nt	Article No.	contact Merca	do for more informatior	۱.
11	Side support cushion X	S, 15x8 cm	276463			
12	Side support cushion S	, 28x8 cm	276462	16 Seat ErgoMed	ic with hygiene fabric.	To be
13	Side support cushion M	1, 24x14 cm	276461	completed witl	n cover in plush.	
14	Side support cushion L	, 30x16 cm	276460			
Sea	t Slow recovery	Slow recovery	Slow recovery	Slow recovery	Slow recovery	
Size	e in width x depth	ErgoMedic	ErgoMedic coxit	ErgoMedic PLUS	ErgoMedic PLUS c	oxit
15	33x35 cm	155466	-	156466	-	
15	36x39 cm	155465	-	156465	-	
15	36x43 cm	155467	-	156467	-	
15	40x39 cm	155464	165464	156464	166464	
15	40x43 cm	155463	165463	156463	166463	
15	40x48 cm	155469	165469	156469	166469	
15	44x43 cm	155462	165462	156462	166462	
15	44x48 cm	155460	165460	156460	166460	
15	48x48 cm	155461	165461	156461	166461	

MEASUREMENTS ON BACKRESTS

ErgoMedic





Size: 27x32 cm

10 Size: 27x22 cm

9

Art.no: 210467

Art.no: 210466

ErgoMedic PLUS



1	Size: 46x43 cm	Art.no: 251469
2	Size: 42x43 cm	Art.no: 251468
3	Size: 35x43 cm	Art.no: 251463
4	Size: 40x20 cm	Art.no: 251460
5	Size: 35x20 cm	Art.no: 251462



6	Size: 30x34 cm	Art.no: 251465
7	Size: 25x34 cm	Art.no: 251464
8	Size: 23x27 cm	Art.no: 253461
9	Size: 23x20 cm	Art.no: 253460



ErgoMedic PLUS Comfort with wedge



- 1 Size: 46x43 cm
- 2 Size: 42x43 cm
- 3 Size: 35x43 cm
- Art.no: 255469 Art.no: 255468 Art.no: 255463

MEASUREMENTS ON SEATS

ErgoMedic





ErgoMedic Slowrecovery



Art.no: 155463

Size: 40x43 cm

5



9 Size: 33x35 cm Art.no: 155465 Art.no: 155466

ErgoMedic Coxit



1	Size: 48x48 cm	Art.no: 160461
2	Size: 44x48 cm	Art.no: 160460
3	Size: 44x43 cm	Art.no: 160462
4	Size: 40x48 cm	Art.no: 160469
5	Size: 40x43 cm	Art.no: 160463
6	Size: 40x39 cm	Art.no: 160464

ErgoMedic Coxit Slowrecovery



ErgoMedic PLUS





ErgoMedic PLUS Slowrecovery

5

Size: 40x43 cm



Art.no: 156463



. 39 cm 6

13 cm

6	Size: 40x39 cm	Art.no: 156464
7	Size: 36x43 cm	Art.no: 156467
8	Size: 36x39 cm	Art.no: 156465
9	Size: 33x35 cm	Art.no: 156466

ErgoMedic PLUS Coxit



1	Size: 48x48 cm	Art.no: 161461
2	Size: 44x48 cm	Art.no: 161460
3	Size: 44x43 cm	Art.no: 161462
4	Size: 40x48 cm	Art.no: 161469
5	Size: 40x43 cm	Art.no: 161463
6	Size: 40x39 cm	Art.no: 161464

ErgoMedic PLUS Coxit Slowrecovery



1	Size: 48x48 cm	Art.no: 166461
2	Size: 44x48 cm	Art.no: 166460
3	Size: 44x43 cm	Art.no: 166462
4	Size: 40x48 cm	Art.no: 166469
5	Size: 40x43 cm	Art.no: 166463
6	Size: 40x39 cm	Art.no : 166464

Wooden seat disc





1	Size: 48x48 cm	Art.no: 150901
2	Size: 44x48 cm	Art.no: 150900
3	Size: 44x43 cm	Art.no: 150902
4	Size: 40x48 cm	Art.no: 150909
5	Size: 40x43 cm	Art.no: 150903
6	Size: 40x39 cm	Art.no: 150904
7	Size: 36x43 cm	Art.no: 150907
8	Size: 36x39 cm	Art.no: 150905
9	Size: 33x35 cm	Art.no: 150906

1	Size: 48x48 cm	Art.no: 150611
2	Size: 44x48 cm	Art.no: 150610
3	Size: 44x43 cm	Art.no: 150612
4	Size: 40x48 cm	Art.no: 150619
5	Size: 40x43 cm	Art.no: 150613
6	Size: 40x39 cm	Art.no: 150614
7	Size: 36x43 cm	Art.no: 150617
8	Size: 36x39 cm	Art.no: 150615
9	Size: 33x35 cm	Art.no: 150616









EXPLANATION OF THE STRUCTURE OF THE ARTICLE NUMBER Article No: 150460. = "ErgoMedic seat, plush grey 44x48 cm"

Digit 1 = Type of upholstered article (1 = seat, 2 = back panel, etc.). Digits 2-3 = Seat system (ErgoMedic, ErgoMedic PLUS, etc.). Digit 4 = Material (plush, skai, etc). Digit 5 = Colour. Digit 6 = Size.

UPHOLSTERY

66

Plush 41 Bordeaux	42 Brown	43 Green	44 Marine blue	45 REAL blue	46 Grey	48 Black		
Skai 31 Red	32 Brown	35 Blue	36 Grey	38 Black				
Patterned fab 71 Red	r ic, Adult 75 Blue	78 Black	Patterned fabri 81 Red	c, Child 85 Blue 89 Yellow				
Alternative covering60 Hygiene fabric96 ESD antistatic								
Covering66 Incontinence cover, grey plastic-coated towel73 Child/Yellow cotton74 Child/Blue cotton								

Coverings are available for all of Mercado Medic's upholstered parts, visit www.mercado.se or contact us for more information.



Mercado Medic AB Tryffelslingan 14 181 57 Lidingö Sweden Phone & FaxPh+ 46 8 555 143 00Fax+ 46 8 555 143 99

Email & Internet Email info@mercado.se Internet www.mercado.se