



# NORA ECO

Effective support in day-to-day care









# Improving mobility and providing a higher quality of life.

The stand-up and raising aid NORA ECO was specially designed for residents who still have the ability to stand for a short amount of time with the help of supporting tools. With NORA ECO caregivers have the ideal tool to transfer these residents safely in a back-friendly and ergonomical way. It also offers the excellent opportunity to mobilize your residents, stabilize and improve their musculoskeletal system.

NORA ECO ist easy to maneuver and can be used by one caregiver. BEKA Hospitec also offers suitable stand-up slings.



#### Improves the mobility of your residents

Residents who are being transferred with the NORA ECO learn to stabilize both the body and help circulation during mobilization. This can lead to a better quality of life for your resident.



## Comfortable, flexible knee pads

The ergonomically shaped knee pads are individually adjustable in height to perfectly suit each resident. They prevent painful pressure on the kneecaps during the standing up motion.



## Easy to maneuver by only one caregiver

The highly maneuverable mobility aid can be conveniently operated by one caregiver. Thanks to the smart ergonomic concept, back injuries and incorrect strain on caregivers can be prevented.

#### Technical Data

Length	1.014 mm	40 inch
Outside width legs, chassis closed	702 mm	27.6 inch
Inside width legs, chassis spread	870 mm	34.2 inch
Weight	approx. 51 kg	112 lb
Height (min. / max.)	1.170 mm / 1.675 mm	46 / 66 inch
SWL	200 kg	440 lb
Battery	24 V	

#### **Product Information**

Electrically, spreadable chassis	
Flexible, height-adjustable knee pads with safety belts	
Emergency lowering function	
Ergonomic handles	
Hand control with symbols for all functions and display of the battery status	

## Available Options



24 V replacement battery



Wall charger

## Sit-to-stand slings



Sit-to-stand sling with clip

Arm and back support with anti-slip coating Double safety straps for the chest area Additionally padded chest strap

#### **BEKA Hospitec GmbH**

Am Rübenmorgen 3 35582 Wetzlar-Dutenhofen info@beka-hospitec.de



- T +49 641 922 20-0
- Fax +49 641 922 20-20
- www.beka-hospitec.com