

Multipurpose patient trolley

Comfort and reliability in acute care and recovery



Emergo patient trolley

The classic Emergo patient trolley is available in a range of models with either 2 or 3-sectional mattress base. The safe and sturdy construction makes it ideal for heavy use. The versatile range of options and adjustments ensures comfort, safety and reliability in various situations such as emergency care and ICU, patient transportation, recovery and minor procedures. Designed and manufactured in Finland, Emergo is a long lasting, good-value patient trolley for any clinic or hospital.

Cost-efficient and long-lasting


- Modular structure, can easily be equipped according to your needs
- A wide range of accessories available
- · Sturdy construction, materials resistant to wear and tear

Comfort and safety for patients

- Comfortable thanks to adjustable mattress base
- Can be lowered to an adequate height for improved patient comfort
- Optional storage compartments to keep the patient's personal items close at hand
- Side rails enhance patient safety
- Also suitable for heavy patients

Ergonomic and convenient for health care staff

- · Light, gas spring-assisted adjustments
- Mattress base rises to an ergonomic working height
- · Side rails are adjustable with a single hand
- Easy to clean, also the X-ray mattress base
- Great to manoeuvre even in tight spaces, fifth wheel improves mobility even further



Factory-assembled options and accessories



Castors 1. Design Ø 150 or 200 mm, antistatic 2. PJP Ø 150 mm

Fixed, without castors (eg. for ships).

Storage



ABS storage box 100006620 large 100006622 large, for dual pedals 100006625 large, for fifth wheel 100006626 small, only with large storage box

(Not compatible with fixed height frame.)



Wire basket under back section 582 100000582


X-ray compatible mattress base, easy to clean


Available in models 6270 and 6280




 $5^{\mbox{th}}$ wheel for improved mobility

Wire basket under foot section 583 100000583





On the left is shown how the turning



Accessory basket 9784 128009784

Accessory holders and rails



Paper roll holder 733 100000733

Infusion rods



Accessory rail (pair) 700 (replaces side rail, arm rest in picture not included) 100000700



Oxygen bottle holder 589 100000589 for 5 l bottles 100005860 for 2 l bottles



Infusion rod, 4-hooks 499 100000499



Fixed, foldable infusion rod for subframe 6678 100006678

Other accessories



Hygienic and antistatic mattresses Size 650mm x 1950mm Thickness 50mm, 80mm or visco-elastic 100mm



Patient restraining straps (set) 790 100000790



Side rails (pair) 100006713 chromed/yellow



Push bar, foldable 6791 100006793 chromed/yellow 100006795 epoxy coated/yellow

Emergo Models

	Emergency Transport Minor procedures	Emergency ICU Minor procedures		Emergency Transport Minor procedures Recovery
Model	6250	6270	6280	6350
Mattress base sections	2	2	2	3
X-ray compatible			\checkmark	
Trendelenburg & Reverse Trendelenburg	~	~	~	~
Adjustable height			\checkmark	
Subframe			\checkmark	
5th wheel	optional	optional	optional	optional

These are suggested trolley models for the quoted applications. Other configurations are also available. Contact your Merivaara representative for further details.

Technical Specifications

Mattress base	2 or 3-sectional, chromed or epoxy-coated steel frame		
Mattress base width	705 mm		
Mattress base length	2075 mm, with subframe 2150 mm		
Recommended mattress size	650 mm x 1950 mm		
Height adjustment min – max, hydraulic	560 mm – 930 mm (with 150 mm castors) 610 mm – 980 mm (with 200 mm castors) X-Ray mattress base increases height 60 mm		
Height in fixed model (100006230)	740 mm (with 150 mm castors)		
Back section angle, gas-spring assisted	0° – 70°		
Leg section angle, gas-spring assisted	0° – 45° (in 3-sectional model)		
Trendelenburg, gas-spring assisted	14°		
Reverse Trendelenburg, gas-spring assisted	7°		
SWL (safe working load)	225 kg		

Custom-made configurations, width and height adjustments are available on special request, contact your nearest Merivaara representative for further details.



Merivaara Corp. Puustellintie 2, FI-15150 Lahti, FINLAND Tel. +358 3 3394 611, fax +358 3 3394 6144 merivaara@merivaara.com www.merivaara.com

C€

The status of the information, specifications and illustrations in this publication is indicated by the date of publication. Whilst the greatest care has been taken over its compilation, Merivaara Corp.is not responsible for any errors or omissions. Merivaara Corp. reserves the right to make changes to the technology, features, specifications and design of the equipment and models without notice. Please check the Merivaara website for the latest updated specifications and other information.

All trademarks are Merivaara property unless otherwise stated. Patents granted or pending.

©Merivaara Corp. DOC ID: 6052EN-3, date 10/2016.

